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Capsular Tension Ring

Contains:

One blue sterile PMMA capsular tension ring.

Description:

Acrivet capsular tension rings are made of ultrapure PMMA.

Capsular tension rings are used after extracapsular cataract extraction in combination with intraocular lens implantation.

Sterilization:

Acrivet capsular tension rings are sterilized with ethylene oxide. A validated degassing process guarantees the elimination of ethylene oxide residues.

Packaging:

The capsular tension ring is packaged sterile and dry in a primary packaging contained in a sterile transparent pouch (Medipeel pouch).

Handling:

To remove the capsular tension ring, open the Medipeel pouch at the defined spot.

Indications:

- Capsular tension rings are used for
- Circular tensioning of the capsular bag
 - Stabilization of the capsular bag in cases of zonular defects, unless the defect is not bigger than 120°
 - Prevention of intraocular lens luxation
 - Stabilization during implantation of foldable lenses

Capsular tension rings can be used in combination with intraocular lenses.

Contraindications:

- Chronic uveitis
- Damaged capsular bag
- Acute inflammation
- Uncontrolled chronic glaucoma

Interactions:

No interactions with drugs are known.

Precautions:

1. Do not resterilize capsular tension rings.
2. Do not implant capsular tension rings from damaged Medipeel pouches.
3. Never use capsular tension rings beyond the expiration date.

4. Store capsular tension rings away from moisture and light. Store capsular tension rings between 2 °C and 50 °C. Do not use capsular tension rings from a moist package.
5. Surgeon should be trained for this kind of operation before implant capsular tension rings.

Return of unsterile or damaged products:

Must be returned in original packaging with capsular tension ring and its serial number.

Contaminated products should not be returned.

Disclaimer of warranty:

Valeant Med Sp. z o.o. is neither responsible for the misuse of the product or in the choice of the capsular tension ring for the respective patient or for use in the appropriate clinical situation of the patient.

Used symbols:

	Sterilized with ethylene oxide		Store away from sunlight
	Do not reuse		Store dry
	Serial number		Observe instruction for use
	Use before		Do not resterilize

BAUSCH + LOMB

 Valeant Med Sp. z o.o.
ul. Ryzowa 31, 02-495 Warszawa, Poland
Tel: +48 22 578 16 00, Fax: +48 22 578 16 18

Acrivet is a trademark of Valeant Pharmaceuticals International, Inc. or its affiliates.
©Bausch & Lomb Incorporated.

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